

**Lactanet - Valacta**

555, boul. des Anciens-Combattants  
Sainte-Anne-de-Bellevue, QC H9X 3R4  
T 514-459-3030 ext 7990  
Email: microlab@lactanet.ca  
1-800-bon-lait



ATTN: Lab manager / supervisors

Date: 2020-09-14

Instrumental counts of canadian fluoro-opto-electronic counters  
Batch : SCC-09-2020

| Standard        | Results (count/ml) |         |         |
|-----------------|--------------------|---------|---------|
|                 | lowest             | average | highest |
| Low (blue)      | 106 000            | 135 000 | 159 000 |
| Medium (yellow) | 468 000            | 516 000 | 563 000 |
| High (red)      | 856 000            | 942 000 | 994 000 |

For best results:

- Expiry date: 2020-10-20
- Do not subdivide samples into smaller portion
- Warm to ± 40°C and shake gently to disperse butterfat
- Discard used vials
- Do not re-heat samples
- Do not use leaking vials

|  |   |  |
|--|---|--|
|  |   |  |
| <input checked="" type="checkbox"/> M-H. Castonguay, M.Sc., Mcb.A. | <input type="checkbox"/> K. M. Amine, M.Sc., Mcb.A. |  |

|   |
|---|
| Cette version remplace et annule toute version antérieure, le cas échéant / If applicable, this version replaces all previous issues.   |
| Ce document est confidentiel / This document is confidential.   |
| Les résultats et commentaires, le cas échéant, ne se rapportent qu'aux échantillons décrits plus haut / Results and comments, if applicable, are only related to the above samples.                         |
| Ce certificat ne doit pas être reproduit, sinon en entier, sans l'autorisation écrite du laboratoire / This certificate shall not be reproduced except in full, without written approval of the laboratory. |

Fin du rapport / End of report